A Patient's Guide to the VenaTech® Convertible Vena Cava Filter

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en CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician

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en **WARRANTY:** B. BRAUN INTERVENTIONAL SYSTEMS INC. WARRANTS THAT THIS PRODUCT WAS MANUFACTURED ACCORDING TO APPLICABLE STANDARDS AND SPECIFICATIONS. PATIENT CONDITION AND CLINICAL TREATMENT MAY AFFECT THE PERFORMANCE OF THIS PRODUCT. USE OF THIS PRODUCT SHOULD BE IN ACCORDANCE WITH THE INSTRUCTIONS PROVIDED AND AS DIRECTED BY THE PRESCRIBING PHYSICIAN.

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Inside This Book	Page
GLOSSARY.....	4
PULMONARY EMBOLISM AND VENA CAVA FILTERS.....	4
PULMONARY EMBOLISM (PE).....	5
WHAT TYPES OF TREATMENT ARE USED FOR PULMONARY EMBOLISM?.....	5
WHAT IS A VENA CAVA FILTER?.....	5
THE VENATECH® CONVERTIBLE™ VENA CAVA FILTER SYSTEM.....	5
INDICATIONS FOR USE.....	6
CONTRAINDICATIONS.....	6
THE IMPLANTATION PROCEDURE.....	7
HOW WILL THE FILTER BE INSERTED?.....	7
HOW LONG DOES THE PROCEDURE USUALLY TAKE?.....	7
WILL I EXPERIENCE DISCOMFORT DURING AND AFTER THE IMPLANT PROCEDURE?.....	7
POTENTIAL ADVERSE EVENTS / TECHNICAL COMPLICATIONS.....	8
AFTER THE PROCEDURE.....	9
MRI SAFETY INFORMATION.....	9
UNDER WHAT CIRCUMSTANCES SHOULD I CONTACT THE DOCTOR RIGHT AWAY?.....	10
RESUMING YOUR NORMAL LIFESTYLE.....	10
CONVERSION PROCEDURE.....	10
QUESTIONS OR CONCERNS.....	11

This brochure is intended to provide you with general information to discuss with your doctor. It is not intended to provide medical care or treatment. You should consult with your doctor regarding the diagnosis or treatment of your medical condition.

GLOSSARY

Anticoagulant - blood thinning medication commonly used to inhibit clots.

Atrial fibrillation - abnormal contractions of the heart muscle resulting in inefficient pumping of blood.

Cardiac arrhythmia - an abnormal heartbeat.

Contraindications - patient attributes that prohibit use of the device.

Deep Vein Thrombosis (DVT) - a blood clot in the lower leg that can partly block the flow of blood.

Embolism - a blood clot that dislodges from one part of the body and travels to another part.

Hematoma - clotted blood located within body tissue.

Hemorrhage - blood leakage from a ruptured blood vessel.

Indications - patient attributes that permit the use of the device.

Inferior Vena Cava (IVC) - the large vein that carries blood from the lower extremities back to the heart and lungs.

Percutaneous - access to the human body through a puncture of the skin.

Pulmonary Embolism (PE) - the condition where a blood clot dislodges and travels through the body and into the lungs, reducing the lung's ability to provide oxygen to your body.

Thrombus - a blood clot.

Thrombosis - a blood clot that forms within a blood vessel or the heart.

Thromboembolism - a blood clot that dislodges from one part of the body and travels to plug another vessel.

Vena Cava Filter - An expandable metal device designed to trap blood clots before they reach the lungs.

PULMONARY EMBOLISM AND VENA CAVA FILTERS

Your doctor has given you this booklet to help you learn more about pulmonary embolism – what causes it, how it can affect your body and, most importantly, how it can be treated. After reading the booklet, talk to your doctor about any questions you have. It is important to remember that each patient is different and that only your doctor can give you information about the details of your specific treatment.

PULMONARY EMBOLISM (PE)

When blood clots form in the lower limbs, they can partly block the flow of blood from the legs. This condition is called deep vein thrombosis or DVT.

If a blood clot dislodges and travels through a person's body to his or her lungs, it causes a condition called pulmonary embolism or PE. A pulmonary embolism reduces or blocks your lung's ability to supply oxygen to your body.

A pulmonary embolism can cause shortness of breath with mild exertion, and can impair lung function. It may even cause heart failure or death. According to the Centers for Disease Control, between 60,000 and 100,000 people in the U.S. die each year from pulmonary embolism.

WHAT TYPES OF TREATMENT ARE USED FOR PULMONARY EMBOLISM?

The most common treatment is a group of medications called anticoagulants or "blood thinners." However, there are some patients who, for a variety of medical reasons, cannot take anticoagulants. For these individuals, a vena cava filter may offer an effective preventative solution.

WHAT IS A VENA CAVA FILTER?

A vena cava filter is an expandable metal device specially designed to trap blood clots before they reach the lungs. The filter is placed in the inferior vena cava (IVC) – the large vein that carries blood from the lower extremities back to the heart and lungs – and remains in place to trap clots before they move further up toward the lungs.

THE VENATECH® CONVERTIBLE™ VENA CAVA FILTER SYSTEM

The VenaTech Convertible Vena Cava Filter System consists of an implantable vena cava filter and implantation accessories. The VenaTech Convertible filter is a flexible, symmetrical, cobalt-chromium metal alloy, self-expanding vena cava filter intended for permanent implantation into the infrarenal inferior vena cava to provide protection against pulmonary embolism. When clinically indicated, the clot-trapping features of the VenaTech Convertible filter can be deactivated by percutaneously converting the filter to an open configuration.

INDICATIONS FOR USE

The VenaTech® Convertible™ Vena Cava Filter System is indicated for the prevention of recurrent pulmonary embolism via placement in the vena cava in the following situations:

- Pulmonary thromboembolism when anticoagulants are contraindicated;
- Failure of anticoagulant therapy in thromboembolic diseases;
- Emergency treatment following massive pulmonary embolism where anticipated benefits of conventional therapy are reduced; and
- Chronic, recurrent pulmonary embolism where anticoagulant therapy has failed or is contraindicated.

When clinically indicated, after implantation the VenaTech® Convertible™ filter can be converted to an open configuration to discontinue filtration according to the 'Conversion Procedure' in the Instructions For Use.


CONTRAINDICATIONS

The VenaTech® Convertible™ Vena Cava Filter System is contraindicated for the following situations:

- Vena cava filters must not be implanted in patients with a risk of septic embolism.
- Patients whose vena cava diameter is greater than (>) 28 mm (due to the risk of device migration).
- Pregnant women or women who are planning to become pregnant

THE IMPLANTATION PROCEDURE

The anatomical sites identified below will provide general guidance on those areas that are important in a vena cava filter implantation procedure.

HOW WILL THE FILTER BE INSERTED?

You will be taken to the operating room or the special procedures laboratory where your physician will choose the appropriate site for access to the vein through which the filter will be placed. The physician passes the filter through a specially designed catheter to the site of deployment in the inferior vena cava, the large central vein that returns blood to the heart.

HOW LONG DOES THE PROCEDURE USUALLY TAKE?

Although it varies depending upon the individual patient and the specific circumstances, the implantation of the filter generally takes less than one (1) hour.

WILL I EXPERIENCE DISCOMFORT DURING AND AFTER THE IMPLANT PROCEDURE?

Local anesthesia, plus a mild sedative that might be taken before the procedure, will normally result in little to no discomfort while the filter is being implanted.

Slight discomfort may be experienced after the procedure, specifically at the site of vessel access.

POTENTIAL ADVERSE EVENTS / TECHNICAL COMPLICATIONS

Possible adverse events include, but are not limited to, the following:

- Access site thrombosis
- Air embolism
- Atrial fibrillation
- Blood flow restriction
- Caval thrombosis or caval occlusion
- Damage to the inferior vena cava
- Death
- Deep vein thrombosis
- Embolization of the device resulting in cardiac arrhythmia or compromise of cardiac valve function
- Extravasation of contrast material at time of vena cavogram causing tissue damage
- Filter fracture
- Filter migration
- Filter tilting
- Hematoma or nerve injury at access or conversion site
- Hemorrhage
- Infection, including access site infection
- Perforation of the inferior vena cava and/or adjacent organs or vertebral bodies
- Pulmonary embolism
- Shortness of breath
- Venous insufficiency

AFTER THE PROCEDURE

How long will the filter last?

The VenaTech Convertible Vena Cava Filter is designed to be a permanent implant and will not need to be removed. When clinically indicated, the clot-trapping features of the VenaTech Convertible filter can be 'converted' to an open configuration by percutaneously unlocking the mechanism that secures the filter head to the filtering legs and removing the filter head.

How long will it take to fully recover?

Recovery from the procedure should be rapid, although the specific length of time will vary from patient to patient, depending upon factors such as age, general state of health, etc.

Can the filter become clogged?

When one or more large, life-threatening blood clots are captured by the filter, the flow of blood through the inferior vena cava can be partially or entirely obstructed. Clots trapped in the filter will normally dissolve as blood flows across them. New or significantly increased swelling in the legs may indicate that the filter has trapped large clots, and if you notice such swelling, you should promptly contact your physician.

MRI SAFETY INFORMATION

Non-clinical testing has demonstrated the VenaTech® Convertible™ Vena Cava Filter is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5T or 3T
- Maximum spatial field gradient of 1,160 gauss/cm (11.6 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

Note: The MRI information is related to the implanted filter only. The implantation accessories were not tested.

Under the scan conditions defined above, the VenaTech Convertible is expected to produce a maximum temperature rise of less than 3.3° C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 11 mm from the VenaTech Convertible when imaged with a gradient echo pulse sequence and a 3T MRI system.

Patient Registration: Patients may consider registering the MR conditions under which the VenaTech® Convertible™ filter can be safely scanned with the MedicalAlert Foundation (888-633-4298 or www.medicalert.org).

UNDER WHAT CIRCUMSTANCES SHOULD I CONTACT THE DOCTOR RIGHT AWAY?

You should contact your physician right away if you experience any of the following:

- Sudden onset of chest pain accompanied by shortness of breath
- Swelling in both legs
- Unexplained pain in the abdomen

RESUMING YOUR NORMAL LIFESTYLE

Should I restrict my activities after the filter implantation or procedure?

The implantation of a vena cava filter is not necessarily a reason to restrict your normal activity level; however, each patient is unique and there may be other medical reasons for doing so. Be sure to discuss with your doctor what level of activity is most appropriate for you following the implantation procedure.

CONVERSION PROCEDURE

In cases where the physician determines that after filter implantation it is clinically indicated to discontinue filtration, the VenaTech Convertible filter can be converted to an open configuration to no longer function as a filter. Filter conversion is achieved by percutaneously unlocking the mechanism that secures the filter head to the filtering legs and removing the filter head.

VenaTech® Convertible™ Vena Cava Filter es VenaTech® Convertible™ Filtro de Vena Cava

- en** CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician
- es** PRECAUCIÓN: La ley federal (EE. UU.) establece que este dispositivo solo se puede vender por o bajo la orden de un médico.
- en** **WARRANTY:** B. BRAUN INTERVENTIONAL SYSTEMS INC. WARRANTS THAT THIS PRODUCT WAS MANUFACTURED ACCORDING TO APPLICABLE STANDARDS AND SPECIFICATIONS. PATIENT CONDITION AND CLINICAL TREATMENT MAY AFFECT THE PERFORMANCE OF THIS PRODUCT. USE OF THIS PRODUCT SHOULD BE IN ACCORDANCE WITH THE INSTRUCTIONS PROVIDED AND AS DIRECTED BY THE PRESCRIBING PHYSICIAN.
- es** **GARANTÍA:** B. BRAUN INTERVENTIONAL SYSTEMS INC. GARANTIZA QUE ESTE PRODUCTO FUE FABRICADO DE ACUERDO CON LAS ESPECIFICACIONES Y NORMATIVAS APLICABLES. EL CUADRO Y EL TRATAMIENTO CLÍNICO DEL PACIENTE PUEDEN AFECTAR EL RENDIMIENTO DE ESTE PRODUCTO. ESTE PRODUCTO DEBE UTILIZARSE DE ACUERDO CON LAS INSTRUCCIONES SUMINISTRADAS Y CON LO INDICADO POR EL MÉDICO.

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B. Braun Interventional Systems Inc.